







# TRIMION AFRICA™

## PICO LASER PROFESSIONAL SYSTEM

Ink is messy; a single wavelength can't tackle everything. That's why this advanced pico laser system is equipped with 1064nm technology for deeper dermal targeting and darker skin tones, 532nm for red, orange and warm pigments, plus optional specialized probes for more challenging pigment tones. One professional platform — a complete spectrum aesthetic solution.

### TREATMENT PROBES



	<b>532NM PROBE</b>	Epidermal Pigment Lesions, Red Yellow, Coffee Tattoo Removal.
	<b>1064NM PROBE</b>	Honeycomb Shape, Break The Melanin More Dispersed!
	<b>532NM-F50</b>	White Porcelain Doll, Remove Facial Color Spots, Make The Skin Feel Clean And Fresh.
	<b>1064NM-F50</b>	Dermal Pigment Lesions And Dark Color Tatoo Removal.
	<b>1320NM-F0</b>	Hollywood Peel Treatment, Carbon Peel Treatment.

### Core Professional Functions

- Tattoo removal support for multi-colour pigment treatments
- Pigmentation-focused treatments including freckles, age spots and sun spots
- Skin rejuvenation and collagen stimulation support
- Nevus & birthmark appearance lightening support
- Carbon peel treatment support for pores and refreshed skin texture

## Treatment Areas & Before / After

### TREATMENT SITES

Suit for different body treated area  
more faster treatment




FACE



ARMPIT



ARM



BIKINI LINE



BACK



LEG

### BEFORE & AFTER



BEFORE

Thick and thick  
armpit hair



AFTER

Thinning of armpit  
hair and gradual  
reduction of dark  
spots



BEFORE

Thick leg hair



AFTER

slowed hair growth



BEFORE

hair is  
clearly visible



AFTER

Clean and smooth,  
get rid of  
embarrassment



BEFORE

clear beard



AFTER

Pores are reduced,  
skin texture is refined

## TREATMENT SITES

Suit for different body treated area  
more faster treatment



FACE



ARMPIT



ARM



BIKINI LINE



BACK




LEG

## BEFORE & AFTER



# Professional Probe Configuration



## Advanced Wavelength Configuration

- 1064nm Probe — Designed for deeper dermal pigment support and darker pigment targeting
- 532nm Probe — Focused on red, orange and lighter pigment shades
- 1320nm Carbon Peel Head — Professional carbon facial treatment support
- Optional probe configurations available for expanded pigment versatility

## Suggested Operational Parameters

Treatment Mode	Frequency	Distance	Energy	Beam	Interval
Eyebrow Removal	1-3Hz	3-5cm	200-400	5-8mm	3-4w
Eye Line Removal	1-3Hz	3-5cm	200-400	5-8mm	3-4w
Lip Line Removal	1-3Hz	3-5cm	300-500	5-8mm	3-4w
Tattoo Removal	1-3Hz	3-5cm	300-500	5-8mm	3w
Spot & Naevi	1-3Hz	3-5cm	300-500	5-8mm	3w
Freckle/Birthmark	1-3Hz	3-5cm	300-500	5-8mm	3w

**TRIMION AFRICA™**

Professional Aesthetic Equipment Supplier

WhatsApp: 061 118 5978

117 Stephen Dlamini Road, Overport, Durban

Professional-use systems only. Results may vary depending on operator protocols and individual response.